

## Product Datasheet

### IL-37 Recombinant Protein (orb1230437)

<b>Description</b>	IL-37 Recombinant Protein
<b>Target</b>	IL37
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	The recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted protein is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing the recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	MVHTSPKVKN LNPKKFSIHD QDHKVLVLD GNLIAVPDKN YIRPEIFFAL ASSLSSASAE KGSPILLGVS KGEFCLYCDK DKGQSHPSLQ LKKEKLMKLA AQKESARRPF IFYRAQVGSW NMLESAAHPG WFICTSCNCN EPVGVTDKFE NRKHIEFSFQ PVCKAEMSPS EVSD
<b>Purity</b>	≥ 95% by SDS-PAGE gel and HPLC analyses.
<b>Source</b>	E. coli
<b>Uniprot ID</b>	<a href="#">Q9NZH6</a>
<b>Expiration Date</b>	6 months from date of receipt.