

Product Datasheet

TNF-alpha Recombinant Protein (orb1230507)

Description	TNF-alpha Recombinant Protein
Target	TNF
Form/Appearance	Lyophilized
Storage	The lyophilized TNF-alpha recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted TNF-alpha is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing TNF-alpha recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Note	For research use only
Protein Sequence	VRSSRTPSD KPVAHVVANP QAEGQLQWLN RRANALLANG VELRDNQLVV PSEGLYLIYS QVLFKGGQCP STHVLLHTI SRIAVSYQTK VNLLSAIKSP CQRETPEGAE AKPWYEPIYL GGVFQLEKGD RLSAEINRPD YLDFAESGQV YFGIIAL
Purity	Greater than 98% by SDS-PAGE gel and HPLC analyses.
Source	E. coli
Uniprot ID	P01375
NCBI	NP_000585.2
Expiration Date	6 months from date of receipt.