

## Product Datasheet

### IFN- $\gamma$ Recombinant Protein (orb1230559)

<b>Description</b>	IFN- $\gamma$ Recombinant Protein
<b>Target</b>	IFNG
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	The lyophilized IFN- $\gamma$ recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted IFN- $\gamma$ stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing IFN- $\gamma$ recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	MQDPYVKEAE NLKKYFNAGH SDVADNGTLF LGILKNWKEE SDRKIMQSQI VSFYFKLFKN FKDDQSIQKS VETIKEDMNV KFFNSNKKKR DDFEKLTNYS VTDLNVQRKA IHELIQVMAE LSPAAKTGKR KRSQMLFQGR RASQ
<b>Purity</b>	Greater than 98% by SDS-PAGE gel and HPLC analyses.
<b>Source</b>	E. coli
<b>Uniprot ID</b>	<a href="#">P01579</a>
<b>NCBI</b>	<a href="#">NP_000610.2</a>
<b>Expiration Date</b>	6 months from date of receipt.