

Product Datasheet

IL-9 Recombinant Protein (orb1230580)

Description	IL-9 Recombinant Protein
Target	IL9
Form/Appearance	Lyophilized
Storage	The lyophilized IL-9 recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted IL-9 is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing IL-9 recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Note	For research use only
Protein Sequence	MQRCSTSWG I QHTSYLIENL KDDPSSK CSC SANVT SCLCL PIPSD DCTTP CFQEGMSQVT NATQQSKFSP FFFRVKRIVE TLKSNKCQFF SCEKPCNQTT AGNTVSFLKS LLKTFQKTEV QVQRSRA
Purity	Greater than 98% by SDS-PAGE gel and HPLC analyses.
Source	E. coli
Uniprot ID	D4A819
NCBI	NP_001099217.1
Expiration Date	6 months from date of receipt.