

Product Datasheet

MMP-3 Recombinant Protein (orb1230591)

Description

MMP-3 Recombinant Protein

Target

MMP3

Form/Appearance

Lyophilized

Storage

The lyophilized MMP-3 recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted MMP-3 is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing MMP-3 recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.

Note

For research use only

Protein Sequence

MRTFPGIPKW RKTHLTyrIV NYTPDLPKDA
VDSAVEKALK VWEEVTPLTF SRLYEGEADI
MISFAVREHG DFYPFDGPGN VLAHAYAPGP
GINGDAHFD DDEQWTKDTTG TNLFLVAAHE
IGHSLGLFHS ANTEALMYPL YHSLDLTRF
RLSQDDINGI QSLYGPPDS PETPLVPTEP
VPPEPGTPAN CDPALSFDV STLGEILIF
KDRHFWRKSL RKLEPELHLI SSFWPSLPSG
VDAAEYVTSK DLVFIFKGNQ FWAIRGNEVR
AGYPRGIHTL GFPPTVRKID AAISDKEKNK
TYFFVEDKYW RFDEKRNSME PGFPKQIAED
FPGIDSKIDA VFEEFGFFYF FTGSSQLEFD
PNAKKVTHTL KSNSWLNC

Purity

Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).

Source

E. coli

Uniprot ID

P08254

NCBI

NP_002413.1

Expiration Date

6 months from date of receipt.