
Product Datasheet

Soluble RANK Ligand Recombinant Protein (orb1230633)

Description

Soluble RANK Ligand Recombinant

| | |
|-------------------------|---|
| Target | Tnfsf11 |
| Form/Appearance | Lyophilized |
| Storage | The lyophilized Soluble RANK Ligand recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted Soluble RANK Ligand is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing Soluble RANK Ligand recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature. |
| Note | For research use only |
| Protein Sequence | PAMMEGSWLD VARRGKPEAQ PFAHLTINAA DIPSGSHKVS LSSWYHDRGW AKISNMTLN GKLRVNQDGF YYLYANICFR HHETSGSVPA DYQLMVYVV KTSIKIPSSH NLMKGGSTKN WSGNSEFHFY SINVGFFFL RAGEEISVQV SNPSLLDPDQ DATYFGAFKV QDID |
| Purity | Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg). |
| Source | E. coli |
| Uniprot ID | Q9ESE2 |
| NCBI | NP_476490.1 |
| Expiration Date | 6 months from date of receipt. |