

Product Datasheet

Soluble RANK Ligand Recombinant Protein (orb1230662)

Description

Soluble RANK Ligand Recombinant

Target	Tnfsf11
Form/Appearance	Lyophilized
Storage	The lyophilized Soluble RANK Ligand recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted Soluble RANK Ligand is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing Soluble RANK Ligand recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Note	For research use only
Protein Sequence	PAMMEGSWLD VAQRGKPEAQ PFAHLTINAA SIPSGSHKVT LSSWYHDRGW AKISNMTLSN GKLRVNQDGF YYLYANICFR HHETSGSVPT DYQLMVYVV KTSIKIPSSH NLMKGGSTKN WSGNSEFHFY SINVGGF FKL RAGEEISIQV SNPSLLDPDQ DATYFGAFKV QDID
Purity	Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
Source	E. coli
Uniprot ID	O35235
NCBI	NP_035743.2
Expiration Date	6 months from date of receipt.