

Product Datasheet

IFN gamma Recombinant Protein (orb1230664)

Description	IFN gamma Recombinant Protein
Target	Ifng
Form/Appearance	Lyophilized
Storage	The lyophilized IFN gamma recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted IFN gamma is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing IFN gamma recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Note	For research use only
Protein Sequence	MHGTVIESLE SLNNYFNSSG IDVEEKSLFL DIWRNWQKDG DMKILQSQII SFYLRLFEVL KDNQAISNNI SVIESHLITT FFSNSKAKKD AFMSIAKFEV NNPQVQRQAF NELIRVVHQL LPESSLRKRK RSRC
Purity	Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
Source	E. coli
Uniprot ID	P01580
NCBI	P01580
Expiration Date	6 months from date of receipt.