

Product Datasheet

TNF alpha Recombinant Protein (orb1230667)

Description	TNF alpha Recombinant Protein
Target	Tnf
Form/Appearance	Lyophilized
Storage	The lyophilized TNF alpha recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted TNF alpha is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing TNF alpha recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Note	For research use only
Protein Sequence	MLRSSSQNSS DKPVAHVVAN HQVEEQLEWL SQRANALLAN GMDLKDNQLV VPADGLYLVY SQVLFGQGQC PDYVLLTHTV SRFAISYQEK VNLLSAVKSP CPKDTPEGAE LKPWYEPIYL GGVFQLEKGD QLSAEVNLPK YLDFAESGQV YFGVIAL
Purity	Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
Source	E. coli
Uniprot ID	P06804
NCBI	NP_001265530.1
Expiration Date	6 months from date of receipt.