

Product Datasheet

IL-3 Recombinant Protein (orb1230681)

Description	IL-3 Recombinant Protein
Target	IL3
Form/Appearance	Lyophilized
Storage	The lyophilized IL-3 recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted IL-3 is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing IL-3 recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Note	For research use only
Protein Sequence	MDTHRLTRTL NCSSIVKEII GKLPEPELKT DDEGPSLRNK SFRRVNLSKF VESQGEVDPE DRYVIKSNLQ KLNCCCLPTSA NDSALPGVFI RDLDDEFKRL RFYMVHLNDL ETVLTSRPPQ PASGSVSPNR GTVEC
Purity	Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
Source	E. coli
Uniprot ID	P01586
NCBI	NP_034686.2
Expiration Date	6 months from date of receipt.