

Product Datasheet

IL-31 Recombinant Protein (orb1230685)

Description	IL-31 Recombinant Protein
Target	IL31
Form/Appearance	Lyophilized
Storage	The lyophilized IL-31 recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted IL-31 is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing IL-31 recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Note	For research use only
Protein Sequence	MTCSLSFGAP ISKEDLRRTI DLLKQESQDL YNNYSIKQAS GMSADESIQL PCFSLDREAL TNISVIIAHL EKVKVLSENT VDTSWVIRWL TNISCFNPLN LNISVPGNTD ESYDCKVFVL TVLKQFSNCM AELQAKDNTT C
Purity	Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
Source	E. coli
Uniprot ID	Q6EAL8
NCBI	NP_083870.1
Expiration Date	6 months from date of receipt.