

Product Datasheet

IL-17F Recombinant Protein (orb1230689)

Description	IL-17F Recombinant Protein
Target	Il17f
Form/Appearance	Lyophilized
Storage	The lyophilized IL-17F recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted IL-17F is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing IL-17F recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Note	For research use only
Protein Sequence	MRKNPKAGVP ALQKAGNCPP LEDNTVRVDI RIFNQNGIS VPREFQNRSS SPWDYNITRD PHRFPSEIAE AQCRHSGCIN AQGQEDSTMN SVAIQQEILV LRREPQGCSN SFRLEKMLLK VGCTCVKPIV HQAA
Purity	Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
Source	E. coli
Uniprot ID	Q7TNI7
NCBI	NP_665855.2
Expiration Date	6 months from date of receipt.