

Product Datasheet

LD78 beta Recombinant Protein (orb1230707)

Description	LD78 beta Recombinant Protein
Target	CCL3L1
Form/Appearance	Lyophilized
Storage	The lyophilized LD78 beta recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted LD78 beta is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing LD78 beta recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Note	For research use only
Protein Sequence	APLAADTPTA CCFSYTSRQI PQNFIADYFE TSSQCSKPSV IFLTKRGRQV CADPSEEWVQ KYVSDLELSA
Purity	Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
Source	E. coli
Uniprot ID	P16619
NCBI	NP_066286.1
Expiration Date	6 months from date of receipt.