

Product Datasheet

GRO gamma Recombinant Protein (orb1230718)

Description	GRO gamma Recombinant Protein
Target	CXCL3
Form/Appearance	Lyophilized
Storage	The lyophilized GRO gamma recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted GRO gamma is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing GRO gamma recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Note	For research use only
Protein Sequence	ASVVTCLRCQ CLQTLQGIHL KNIQSVNVRS PGPHCAQTEV IATLKNGKKA CLNPASPMVQ KIIEKILNKG STN
Purity	Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
Source	E. coli
Uniprot ID	P19876
NCBI	NP_002081.2
Expiration Date	6 months from date of receipt.