

Product Datasheet

RANTES Recombinant Protein (orb1230747)

Description	RANTES Recombinant Protein
Target	CCL5
Form/Appearance	Lyophilized
Storage	The lyophilized RANTES recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted RANTES is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing RANTES recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Note	For research use only
Protein Sequence	SPYSSDTTPC CFAYIARPLP RAHIKEYFYT SGKCSNPAVV FVTRKNRQVC ANPEKKWVRE YINSLEMS
Purity	Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
Source	E. coli
Uniprot ID	P13501
NCBI	NP_002976.2
Expiration Date	6 months from date of receipt.