

Product Datasheet

BLC Recombinant Protein (orb1230759)

| | |
|-------------------------|---|
| Description | BLC Recombinant Protein |
| Target | Cxcl13 |
| Form/Appearance | Lyophilized |
| Storage | The lyophilized BLC recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted BLC is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing BLC recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature. |
| Note | For research use only |
| Protein Sequence | ILEAHYTNLK CRCSGVISTV VGLNIIDRIQ VTPPGNGCPK TEVVIWTKMK KVICVNPRAK WLQRLLRHVQ SKLSSTPQA PVSKRRAA |
| Purity | Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg). |
| Source | E. coli |
| Uniprot ID | O55038 |
| NCBI | NP_061354.1 |
| Expiration Date | 6 months from date of receipt. |