

Product Datasheet

IL-17A Recombinant Protein (orb1230843)

Description	IL-17A Recombinant Protein
Target	IL17A
Form/Appearance	Lyophilized
Storage	The lyophilized IL-17A recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted IL-17A is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing IL-17A recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Note	For research use only
Protein Sequence	MIVKAGITIP RNP GCPNSED KNFPRTVMVN LNIHNRNTNT NPKRSSDYNN RSTSPWNLHR NEDPERYPSV IWEAKCRHLG CINADGNVDY HMNSVPIQQE ILVLRREPPH CPNSFRLEKI LVS VGCTCVT PIVHHVA
Purity	Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
Source	E. coli
Uniprot ID	Q16552
NCBI	NP_002181.1
Expiration Date	6 months from date of receipt.