

Product Datasheet

IL-13 Variant Recombinant Protein (orb1230847)

Description	IL-13 Variant Recombinant Protein
Target	IL13
Form/Appearance	Lyophilized
Storage	The lyophilized IL-13 Variant recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted IL-13 Variant is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing IL-13 Variant recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Note	For research use only
Protein Sequence	MSPGPVPPST ALRELIEELV NITQNQKAPL CNGSMVWSIN LTAGMYCAAL ESLINVSGCS AIEKTQRMLS GFCPHKVSAG QFSSLHVRDT KIEVAQFVKD LLLHLKKLFR EGQFN
Purity	Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
Source	E. coli
Uniprot ID	P35225
NCBI	NP_002179.2
Expiration Date	6 months from date of receipt.