

Product Datasheet

IL-1 Receptor Agonist Recombinant Protein (orb1230860)

Description

IL-1 Receptor Agonist Recombinant

Target	IL1RN
Form/Appearance	Lyophilized
Storage	The lyophilized IL-1 Receptor Agonist recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted IL-1 Receptor Agonist is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing IL-1 Receptor Agonist recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Note	For research use only
Protein Sequence	MRPSGRKSSK MQAFRIWDVN QKTFYLRNNQ LVAGYLQGPN VNLEEKIDVV PIEPHALFLG IHGGKMCLSC VKSGDETRLQ LEAVNITDLS ENRKQDKRFA FIRSDSGPTT SFESAACPGW FLCTAMEADQ PVSLTNMPDE GVMVTKFYFQ EDE
Purity	Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
Source	E. coli
Uniprot ID	P18510
NCBI	NP_000568.1
Expiration Date	6 months from date of receipt.