

Product Datasheet

sDLL-4 Recombinant Protein (orb1230867)

Description

sDLL-4 Recombinant Protein

Target

DLL4

Form/Appearance

Lyophilized

Storage

The lyophilized sDLL-4 recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted sDLL-4 is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing sDLL-4 recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.

Note

For research use only

Protein Sequence

SGVFQLQLQE FINERGVLAS GRPCEPGCRT
FFRVCLKHFQ AVVSPGPCTF GTVSTPVLGT
NSFAVRDDSS GGGRNPLQLP FNFTWPGTFS
LIIEAWHAPG DDLRPEALPP DALISKIAIQ
GSLAVGQNW LDEQSTLTLR LRYSYRVICS
DNYYGDNCSR LCKKRNDHFG HYVCQPDGNL
SCLPGWTGEY CQQPICLSGC HEQNGYCSKP
AECLCRPGWQ GRLCNECIPH NGCRHGTCT
PWQCTCDEGW GGLFCDQDLN
YCTHHSPCKN GATCSNSGQR SYTCTCRPGY
TGVDCELELS ECDSNPCRNG GSCKDQEDGY
HCLCPPGYG LHCEHSTLSC ADSPCFNGGS
CRERNQGANY ACECPPNFTG SNCEKKVDRC
TSNPCANGGQ CLNRGPSRMC RCRPGFTGT
CELVSDCAR NPCAHGGTCH DLENGLMCTC
PAGFSGRRCE VRTSIDACAS SPCFNRATCY
TDLSTDTFVC NCPYGFVGSR CEFPVGLP

Purity

Greater than 95% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).

Source

HEK293 cells

Uniprot ID
[Q9NR61](#)
NCBI
[NP_061947.1](#)
Expiration Date

6 months from date of receipt.