

Product Datasheet

IL-6 Recombinant Protein (orb1230976)

Description	IL-6 Recombinant Protein
Target	Il6
Form/Appearance	Lyophilized
Storage	The lyophilized IL-6 recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted IL-6 is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing IL-6 recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.
Note	For research use only
Protein Sequence	MFPTSQVRRG DFTEDTTPNR PVYTTSQVGG LITHVLWEIV EMRKELCNGN SDCMNDDAL AENNLKLEI QRNDGCYQTG YNQEICLLKI SSGLLEYHSY LEYMKNLKD NKKDKARVLQ RDTETLIHIF NQEVKDLHKI VLPTPISNAL LTDKLESQKE WLRTKTIQFI LKSLEEFKLV TLRSTRQT
Purity	Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
Source	E. coli
Uniprot ID	P08505
NCBI	NP_112445.1
Expiration Date	6 months from date of receipt.