

---

## Product Datasheet

### Visfatin Recombinant Protein (orb1230977)

**Description**

Visfatin Recombinant Protein

**Target**

NAMPT

**Form/Appearance**

Lyophilized

**Storage**

The lyophilized Visfatin recombinant protein is stable for at least 2 years from date of receipt at -20°C. Reconstituted Visfatin is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. As with any protein, exposing Visfatin recombinant protein to repeated freeze / thaw cycles is not recommended. When working with proteins care should be taken to keep recombinant protein at a cool and stable temperature.

**Note**

For research use only

**Protein Sequence**

MPPNTSKVYS YFECREKKTE NSKLRKVKYE  
 ETVFYGLQYI LNKYLKGGVV TKEKIQEAKD  
 VYKEHFQDDV FNEKGWNYIL EKYDGHLP  
 IE IKAVPEGFVI PRGNVLFVTE NTDPECYWLT  
 NWIETILVQS WYPITVATNS REQKILAKY  
 LLETSGNLDG LEYKLHDFGY RGVSSQETAG  
 IGASAHLVNF KGTDTVAGLA LIKKYYGTKD  
 PVPGYSPAA EHSTITAWGK DHEKDAFEHI  
 VTQFSSVPVS VVSDSYDIYN ACEKIWGEDL  
 RHLIVSRSTQ APLIIRPDSG NPLDVLKVL  
 EILGKKFPVT ENSKGYKLLP PYLRVIQGDG  
 VDINTLQEIV EGMKQKMWSI ENIAFGSGGG  
 LLQKLTRDLL NCSFKCSYVV TNGLGINVFK  
 DPVADPNKRS KKGRLSLHRT PAGNFVTL  
 EE GKGDL E EYGQ DLLHTVFKNG KVT  
 KSYSFDE IRKNAQLNIE LEAAHH

**Purity**

Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).

**Source**

E. coli

**Uniprot ID**
[P43490](#)
**NCBI**
[NP\\_005737.1](#)
**Expiration Date**

6 months from date of receipt.