

## Product Datasheet

# SARS-CoV-2 (COVID-19) Spike E Mosaic Recombinant protein (orb1231006)

<b>Description</b>	SARS-CoV-2 (COVID-19) Spike E Mosaic Recombinant protein
<b>Preservatives</b>	CoV 2019 Spike E-Mosaic Protein 1mg/ml solution is supplied in 1x PBS.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Store in working aliquots at -20°C. Avoid freeze/thaw cycles.
<b>Tag</b>	His Tag
<b>Note</b>	For research use only
<b>Protein Sequence</b>	The E.Coli derived recombinant protein contains the Coronavirus 2019 spike Envelope Mosaic immunodominant regions, fused to 6xHis tag at C-terminal.
<b>Purity</b>	> 90% pure as determined SDS-PAGE.
<b>Source</b>	E. coli
<b>Expiration Date</b>	6 months from date of receipt.