


Product Datasheet

GITR / TNFRSF18 Recombinant Protein (orb1231200)

Description

GITR / TNFRSF18 Recombinant Protein



Tested

ELISA, WB

Applications

Target

TNFRSF18

Preservatives

Tris with Glycine, Arginine and NaCl, pH 7.5

Form/Appearance

Lyophilized

Storage

Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.

Tag


Fc Tag

Note


For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 41.0 kDa. The protein migrates as 45-47 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.



Human
GITR, Fc
Tag on SDS-
PAGE under
red...



Immobilized
Human
GITR, Fc
Tag at 2
ug/m...

Protein Sequence

Gln 26 - Glu 161

Purity

> 95% as determined by SDS-PAGE.

Source

HEK293 cells

MW

41.0 kDa

Uniprot ID

[Q9Y5U5](https://www.uniprot.org/entry/Q9Y5U5)

NCBI

[NP_004186.1](https://www.ncbi.nlm.nih.gov/nuccore/NP_004186.1)

Dilution Range

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 41.0 kDa. The protein migrates as 45-47 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Expiration Date

6 months from date of receipt.