


Product Datasheet

TSLP (R127A, R130A) Recombinant Protein (orb1231206)

Description

TSLP (R127A, R130A) Recombinant Protein



Tested

ELISA, WB

Applications

Target

TSLP

Preservatives

PBS, pH 7.4

Form/Appearance

Lyophilized

Storage

Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.

Tag

His Tag


Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 16.7 kDa. The protein migrates as 22-27 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Human
TSLP
(R127A,
R130A), His
Tag on
SD...



Protein Sequence

Tyr 29 - Gln 159

Purity

> 90% as determined by SDS-PAGE.

Source

HEK293 cells

MW

16.7 kDa

Uniprot ID

[Q969D9](#)


NCBI

[R130](#), [NP_149024.1](#), [R127](#)

Dilution Range

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 16.7 kDa. The protein migrates as 22-27 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Immobilized
Human
TSLP R, Fc
Tag at 5
ug...



Expiration Date

6 months from date of receipt.

Immobilized
Human
TSLP
(R127A,
R130A), H...