
Product Datasheet

ICOS / CD278 (C137S, C138S) Recombinant Protein (orb1231356)

Description

ICOS / CD278 (C137S, C138S) Recombinant


Tested Applications

ELISA, WB

Target

Icos

Preservatives

Tris with Glycine, Arginine and NaCl, pH 7.5

Form/Appearance

Lyophilized

Storage

Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.

Tag

Fc Tag

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 40.6 kDa. The protein migrates as 45-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

Glu 21 - Leu 144

Purity

> 95% as determined by SDS-PAGE.

Source

HEK293 cells

MW


40.6 kDa

Uniprot ID
[Q9WVS0](#)
NCBI
[C137](#), [NP_059508.2](#), [C138](#)
Dilution Range


This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 40.6 kDa. The protein migrates as 45-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Expiration Date

6 months from date of receipt.



Mouse ICOS (C137S, C138S), Fc Tag on SDS...



Immobilized Human B7-H2, His Tag at 1 ug...