

Product Datasheet

TSLP (R127A, R130S) Recombinant Protein (orb1231403)

Description

TSLP (R127A, R130S) Recombinant Protein

kDa	M	R
116.0		

Tested Applications

ELISA, WB

Target

TSLP

Preservatives

PBS, pH 7.4

Form/Appearance

Lyophilized

Storage

Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.

Tag

His Tag

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 16.9 kDa. The protein migrates as 23-30 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

Tyr 29 - Gln 159

Purity

> 90% as determined by SDS-PAGE.

Source

HEK293 cells

MW

16.9 kDa

NCBI


[XP_005557555.1](#)

Dilution Range


This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 16.9 kDa. The protein migrates as 23-30 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Expiration Date


6 months from date of receipt.



Cynomolgus
TSLP
(R127A,
R130S), His
Tag ...



Immobilized
Human
TSLP R, Fc
Tag at 5
ug...



Immobilized
Cynomolgus
TSLP
(R127A,
R130S)...