


Product Datasheet

Rhesus macaque Properdin Recombinant Protein (orb1231423)

Description

Rhesus macaque Properdin Recombinant Protein



Tested

WB

Applications

Target

CFP

Preservatives

PBS, pH 7.4

Form/Appearance

Lyophilized

Storage

Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.

Rhesus
macaque
Properdin,
His Tag on
SDS...

Tag

His Tag

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 50.6 kDa. The protein migrates as 55-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

Asp 28 - Leu 469

Purity

> 95% as determined by SDS-PAGE.

Source

HEK293 cells

MW

50.6 kDa

NCBI

[XP_014982784.1](#)

Dilution Range

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 50.6 kDa. The protein migrates as 55-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Expiration Date

6 months from date of receipt.