



Product Datasheet

Streptavidin Recombinant Protein (orb1231439)

Description	Streptavidin Recombinant Protein
Tested Applications	ELISA, WB
Preservatives	NaCl, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	This protein .
Source	E. coli
Dilution Range	This protein .
Expiration Date	6 months from date of receipt.



Immobilized Biotinylated Human CD19, His...



Immobilized Biotinylated SARS-CoV-2 S pr...