


Product Datasheet

HCoV-229E S1 Recombinant Protein (orb1231456)

Description	HCoV-229E S1 Recombinant Protein
Tested Applications	WB
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 58.8 kDa. The protein migrates as 70-110 kDa under reducing (R) condition (SDS-PAGE) due to Glycosylation.
Protein Sequence	Cys 16 - Ala 536
Purity	> 90% as determined by SDS-PAGE.
Source	HEK293 cells
MW	58.8 kDa
Dilution Range	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 58.8 kDa. The protein migrates as 70-110 kDa under reducing (R) condition (SDS-PAGE) due to Glycosylation.
Expiration Date	6 months from date of receipt.



HCoV-229E S1 protein, His Tag on SDS-PAG...