


Product Datasheet

IL-31 RA Recombinant Protein (orb1231526)

Description	IL-31 RA Recombinant Protein
Tested Applications	WB
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized
Storage	Lyophilized Protein should be stored at -20°C or lower for long term storage. Upon reconstitution, working aliquots should be stored at -20°C or -70°C. Avoid repeated freeze-thaw cycles.
Tag	His Tag
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 59.0 kDa. The protein migrates as 90-115 KDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Ala 20 - Glu 519
Purity	> 95% as determined by SDS-PAGE.
Source	HEK293 cells
MW	59.0 kDa
Dilution Range	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 59.0 kDa. The protein migrates as 90-115 KDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Expiration Date	6 months from date of receipt.



Human IL-31 RA, His Tag on SDS-PAGE unde...