


## Product Datasheet

# SARS-CoV-2 (COVID-19) E Recombinant Protein (orb1231545)

<b>Description</b>	SARS-CoV-2 (COVID-19) E Recombinant Protein
<b>Preservatives</b>	Supplied as lyophilized from PBS, pH 7.5. Please reconstitute the protein with ddH2O to the concentration with 1mg/ml.
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt.
<b>Tag</b>	N-His-sumo Tag
<b>Note</b>	For research use only
<b>Application notes</b>	Immunogen
<b>Protein Sequence</b>	Met1-Val75
<b>Purity</b>	> 90% as determined by SDS-PAGE
<b>Source</b>	E. coli
<b>MW</b>	21.04kDa
<b>NCBI</b>	<a href="#">QHD43418.1</a>
<b>Dilution Range</b>	Immunogen
<b>Expiration Date</b>	6 months from date of receipt.



SDS-PAGE  
Image