

## Product Datasheet

### Human ACE2 Recombinant Protein (orb1231558)

<b>Description</b>	Human ACE2 Recombinant Protein
<b>Preservatives</b>	Supplied as lyophilized from TBS, pH 8.0. Please reconstitute the protein with ddH <sub>2</sub> O to the concentration with 1mg/ml.
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week. Store at -20 to -80 °C for twelve months from the date of receipt.
<b>Tag</b>	Without Tag
<b>Note</b>	For research use only
<b>Application notes</b>	Immunogen
<b>Protein Sequence</b>	Met1-Ser740
<b>Purity</b>	> 90% as determined by SDS-PAGE
<b>Source</b>	Mammalian Cells
<b>MW</b>	85.2kDa
<b>NCBI</b>	<a href="#">NP_001358344.1</a> , <a href="#">Q9BYF1</a>
<b>Dilution Range</b>	Immunogen
<b>Expiration Date</b>	6 months from date of receipt.

