

Product Datasheet

Human ACE2 Recombinant Protein (orb1231569)

Description	Human ACE2 Recombinant Protein	Image not found or type unknown
Preservatives	Supplied as lyophilized from PBS, pH 7.5. Please reconstitute the protein with ddH ₂ O to the concentration with 1mg/ml.	SDS-PAGE Image
Form/Appearance	Lyophilized	
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week. Store at -20 to -80 °C for twelve months from the date of receipt.	
Tag	C-His Tag	
Note	For research use only	
Application notes	Immunogen	
Protein Sequence	Ser19-Asp615	
Purity	> 90% as determined by SDS-PAGE quantitative densitometry by Coomassie Blue Staining.	
Source	Mammalian Cells	
MW	71.54kDa	
NCBI	NP_001358344.1 , Q9BYF1	
Dilution Range	Immunogen	
Expiration Date	6 months from date of receipt.	