

Product Datasheet

SARS-CoV-2 (COVID-19) NSP8 Recombinant Protein (orb1231575)

Description SARS-CoV-2 (COVID-19) NSP8 Recombinant Protein

Preservatives Supplied as lyophilized from PBS, pH 7.5. Please reconstitute the protein with ddH₂O to the concentration with 1mg/ml.

Form/Appearance Lyophilized

Storage Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week. Store at -20 to -80 °C for twelve months from the date of receipt.

Tag N-His Tag

Note For research use only

Application notes Immunogen

Protein Sequence Ala1-Gln198

Purity > 85% as determined by SDS-PAGE quantitative densitometry by Coomassie Blue Staining.


Source E. coli

MW 24.18kDa

NCBI [YP_009725304.1](#)

Dilution Range Immunogen

Expiration Date 6 months from date of receipt.



SDS-PAGE
Image