

Product Datasheet

SARS-CoV-2 (COVID-19) NSP10 Recombinant Protein (orb1231581)

| | |
|--------------------------|---|
| Description | SARS-CoV-2 (COVID-19) NSP10 Recombinant Protein |
| Preservatives | Supplied as lyophilized from PBS, pH 7.5. Please reconstitute the protein with ddH2O to the concentration with 1mg/ml. |
| Form/Appearance | Lyophilized |
| Storage | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week. Store at -20 to -80 °C for twelve months from the date of receipt. |
| Tag | N-His Tag |
| Note | For research use only |
| Application notes | Immunogen |
| Protein Sequence | Ala1-Gln139 |
| Purity | > 90% as determined by SDS-PAGE quantitative densitometry by Coomassie Blue Staining. |
| Source | E. coli |
| MW | 17.99kDa |
| NCBI | YP_009725306.1 |
| Dilution Range | Immunogen |
| Expiration Date | 6 months from date of receipt. |



SDS-PAGE Image