

Product Datasheet

SARS-CoV-2 (COVID-19) NSP3 Recombinant Protein (orb1231582)

Description	SARS-CoV-2 (COVID-19) NSP3 Recombinant Protein
Preservatives	Supplied as lyophilized from PBS, pH 7.5. Please reconstitute the protein with ddH ₂ O to the concentration with 1mg/ml.
Form/Appearance	Lyophilized
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week. Store at -20 to -80 °C for twelve months from the date of receipt.
Tag	N-His Tag
Note	For research use only
Application notes	Immunogen
Protein Sequence	Glu1024-Gln1198
Purity	> 90% as determined by SDS-PAGE quantitative densitometry by Coomassie Blue Staining.
Source	E. coli
MW	22.03kDa
NCBI	QHD43415.1
Dilution Range	Immunogen
Expiration Date	6 months from date of receipt.



SDS-PAGE Image