



Product Datasheet

SARS-CoV-2 (COVID-19) Nucleocapsid Recombinant Protein (orb1231585)

Description

SARS-CoV-2 (COVID-19) Nucleocapsid


 SDS-PAGE
Validation
with SARS-
CoV-2
(COV...

 ELISA
Validation
of SARS-
CoV-2
(COVID-
19...

 Western
Blot
Validation
with SARS-
CoV-2 ...

Tested Applications	ELISA, WB
Preservatives	PBS
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	The protein is stable in a liquid state at -70°C for 12 months. Avoid repeated freeze-thaw cycles.
Tag	His tag at C-terminus
Note	For research use only
Protein Sequence	Met1-Ala419
Purity	95%
Source	E. coli
MW	50 kD
Expiration Date	6 months from date of receipt.