

Product Datasheet

Prol1 antibody (orb12316)

Description	Goat polyclonal to Prol1.
Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Mucin 10 / Prol1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments gave an approx 48kDa band in Mouse Eye lysates after 2µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 33.6kDa according to NP_
Protein Sequence	PKRQIISDKTKP
Clonality	Polyclonal
MW	33.6
Dilution Range	ELISA: 1:8000, WB: 2 µg/ml
Expiration Date	12 months from date of receipt.