

Product Datasheet

ACE-2 Recombinant Protein (orb1231623)

Description	ACE-2 Recombinant Protein
Target	ACE2
Preservatives	20mM Tris, 300mM NaCl, 1mM ZnCl ₂ , 10% Glycerol, pH 7.4.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store in working aliquots at -20°C. Avoid freeze/thaw cycles.
Tag	His Tag
Note	For research use only
Protein Sequence	Gln18-Ser740
Purity	> 95% as determined by reducing SDS-PAGE.
Source	HEK293 cells
MW	Mol Mass: 84.63kD; AP Mol Mass: 103kD
NCBI	Q9BYF1
Expiration Date	6 months from date of receipt.