


Product Datasheet

SARS-CoV-2 (COVID-19) ORF8 Recombinant Protein (orb1231626)

Description	SARS-CoV-2 (COVID-19) ORF8 Recombinant Protein
Tested Applications	ELISA, WB
Target	ORF8
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store in working aliquots at -20°C. Avoid freeze/thaw cycles.
Tag	His Tag
Note	For research use only
Protein Sequence	16-121 aa
Purity	> 80% as estimated by SDS-PAGE gel
Source	E. coli
MW	13.4 kD
NCBI	QHD43422.1
Dilution Range	SDS-PAGE
Expiration Date	6 months from date of receipt.



SARS-CoV-2 (COVID-19, 2019-nCoV) ORF8 Re...