

Product Datasheet

PEX10 (CT) Peptide (orb1231632)

Description	PEX10 (CT) Peptide
Tested Applications	Blocking
Target	PEX10
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store PEX10°CT) Peptide at -20°C, stable for one year.
Note	For research use only
Application notes	PEX10 (CT) Peptide is used for blocking the activity of the PEX10 (CT) Antibody
Uniprot ID	O60683
NCBI	O60683
Expiration Date	6 months from date of receipt.