

Product Datasheet

VTI1a Peptide (orb1231755)

Description	VTI1a Peptide
Tested Applications	Blocking
Target	VTI1a
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store VTI1a peptide at -20°C, stable for one year.
Note	For research use only
Application notes	VTI1a peptide is used for blocking the activity of VTI1a antibody.
Uniprot ID	Q96AJ9
NCBI	NP_660207
Expiration Date	6 months from date of receipt.