


Product Datasheet

Nsg1 antibody (orb12318)

Description	Goat polyclonal to Nsg1.
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Target	Nsg1 / Neep21
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Protein Sequence	EKGTKQPLLEDGFD
Clonality	Polyclonal
MW	20.9
Entrez	27065
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml, IHC-P: 1-5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Brain lys...