


Product Datasheet

USP25 Peptide (orb1231923)

Description	USP25 Peptide
Tested Applications	Blocking
Target	USP25
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store USP25 Peptide at -20°C, stable for one year.
Note	For research use only
Application notes	The peptide is used for blocking the activity of USP25 Antibody (catalog number 7765). It usually blocks the antibody activity completely by incubating the peptide with equal volume of antibody for 30 min at 37°C.
Uniprot ID	Q9UHP3
NCBI	NP_037528
Expiration Date	6 months from date of receipt.



Western blot analysis of USP25 in mouse ...