

Product Datasheet

ME2 Peptide (orb1231938)

Description	ME2 Peptide
Tested Applications	Blocking
Target	ME2
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store ME2 peptide at -20°C, stable for one year.
Note	For research use only
Application notes	The peptide is used for blocking the activity of ME2 antibody (catalog number 7723). It usually blocks the antibody activity completely by incubating the peptide with equal volume of antibody for 30 min at 37°C.
Uniprot ID	P23368
NCBI	NP_002387
Expiration Date	6 months from date of receipt.