

Product Datasheet

ME1 Peptide (orb1231939)

Description	ME1 Peptide
Tested Applications	Blocking
Target	ME1
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store ME1 peptide at -20°C, stable for one year.
Note	For research use only
Application notes	The peptide is used for blocking the activity of ME1 antibody (catalog number 7721). It usually blocks the antibody activity completely by incubating the peptide with equal volume of antibody for 30 min at 37°C.
Uniprot ID	P48163
NCBI	NP_002386
Expiration Date	6 months from date of receipt.