

## Product Datasheet

### IL-1F10 Peptide (orb1232004)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | IL-1F10 Peptide   |
| <b>Tested Applications</b> | Blocking  |
| <b>Target</b>              | IL1F10  |
| <b>Preservatives</b>       | PBS pH 7.2 (10 mM NaH <sub>2</sub> PO <sub>4</sub> , 10 mM Na <sub>2</sub> HPO <sub>4</sub> , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide  |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Concentration</b>       | 200 ug/mL   |
| <b>Storage</b>             | Store IL-1F10 peptide at -20°C, stable for one year.  |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | The peptide is used for blocking the activity of IL-1F10 antibody (catalog number 7503). It usually blocks the antibody activity completely by incubating the peptide with equal volume of antibody for 30 min at 37°C. |
| <b>Uniprot ID</b>          | <a href="#">Q8WWZ1</a>  |
| <b>NCBI</b>                | <a href="#">NP_775184</a>   |
| <b>Expiration Date</b>     | 6 months from date of receipt.  |