

Product Datasheet

RBBP8 Peptide (orb1232579)

Description	RBBP8 Peptide
Tested Applications	Blocking
Target	RBBP8
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store RBBP8 peptide at -20°C, stable for one year.
Note	For research use only
Application notes	RBBP8 peptide is used for blocking the activity of RBBP8 antibody.
Uniprot ID	Q99708
NCBI	NP_976037
Expiration Date	6 months from date of receipt.