

Product Datasheet

SYNGR2 Peptide (orb1232665)

Description	SYNGR2 Peptide
Tested Applications	Blocking
Target	SYNGR2
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Form/Appearance	Liquid
Concentration	200 ug/mL
Storage	Store SYNGR2 peptide at -20°C, stable for one year.
Note	For research use only
Application notes	SYNGR2 peptide is used for blocking the activity of SYNGR2 antibody.
Uniprot ID	O43760
NCBI	AAI05993
Expiration Date	6 months from date of receipt.